

- Measure-Up™ Arm-Sleeve
- Silhouette™ Arm Undersleeve


(800)CIRCAID (247-2243)  
 website: www.circaid.com  
 Fax: (858)576-3555

**Color Option for Measure-Up Arm-Sleeve Only**  
 If left blank, foam liner will default to black exterior, beige interior. Shell available in black only.

<b>Date</b>	<b>Reorder</b> <input type="checkbox"/> Yes <input type="checkbox"/> No	<b>Foam Liner Exterior</b> <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Beige <input type="checkbox"/> Turquoise	<b>Foam Liner Interior</b> <input type="checkbox"/> Red <input type="checkbox"/> Beige <input type="checkbox"/> Turquoise
<b>Facility Name</b>		<b>Contact</b>	<b>Patient Name</b>
<b>Billing and/or Shipping Address</b>		<b>Phone Number</b>	
<b>Billing and/or Shipping Address</b>		<b>Fax Number</b>	
<b>City, State, Zip</b>		<b>Account Number</b>	<b>Purchase Order Number</b>

**Shipping Method:**  Ground  3 Day  2 Day  Next Day ( ) Early a.m. ( ) 10:30a.m.

<b>Length</b>	
Wrist to Elbow Crease	
Left	Right



<b>Length*</b>	
Wrist to Axilla	
Left	Right

<b>Circumference Measurements</b>	
centimeters only	
Palm of Hand	
Wrist	
5cm	
10cm	
15cm	
20cm	
25cm	
30cm	
35cm	
40cm	
45cm	
50cm	
55cm	

**Measure in centimeters:**

1. Palm circumference
2. Arm circumference every 5cm starting at wrist
3. Medial length with arm bent at 35° from wrist to axilla
4. Medial length from wrist to elbow crease

**Special Instructions**

- Arm should be at its smallest
  - Maintain arm size at time of measuring to assure a proper fit of garment
  - Facility and patient acknowledge custom-made products are warranted to be free of defects in materials and workmanship and can be returned for repair but not a refund.
- Initials** \_\_\_\_\_

**Finished Length Rules**

Silhouette undersleeve length (measured medially) is determined from the wrist-to-axilla length provided +/- 1cm.

Measure-Up length (measured medially) is determined from wrist-to-axilla length provided. **Note:** If the length is equal to the borderline, sleeve is built to the longer length unless specifically requested to build to the shorter length:

Armsleeve medial length from wrist to top	Min	Max	Number bands on armsleeve
37	33	39.4	6
43	39.4	45.8	7
49	45.8	52.1	8
56	52.1	58.5	9