




Juxta-Lite Legging


The Perfect Alternative to Compression Stockings
Standard Line Measuring and Sizing Chart

Circumference Measurements


Position	
c^C	Calf - Circumference at largest dimension
c^B	Ankle - Circumference at smallest dimension



c^{A1}
Arch Circumference For Juxta-Lite AFW Only



/D
Medial length from floor to the knee crease
Legging Length*
Short : 28cm
Long : 33cm



Patient:

- A. Should be standing.
- B. Limb should be at its smallest and limb size at time of measurements should be maintained to assure a proper fit of garment.
- C. Be sure to select a length (**/D**), legging size (**c^B**, **c^C**), and foot size (**c^{A1}**) if ordering a detached Juxta-Lite™ Ankle-Foot Wrap.

Note: The leggings are also available with a detached one size fits all compression anklet.

Sizing Rules

Measure accurately. If borderline between sizes, consider amount of edema present & anticipated reduction. If significant edema is present and reduction is anticipated, select smaller size; if edema is minimal, select larger size. If measurements do **NOT** fall within **any [one] size category**, a different product type (Juxta-Fit™) will need to be ordered.

	Measurements		Small	Medium	Large	Large Full Calf	XL	XL Full Calf	XXL
	Left	Right							
c^B Ankle			19 - 22cm	22 - 27cm	25 - 32cm	25 - 32cm	30 - 37cm	30 - 37cm	35 - 42cm
c^C Calf			25 - 36cm	33 - 44cm	41 - 51cm	51 - 61cm	37 - 48cm	48 - 61cm	53 - 64cm
/D *Length			Under 44cm Over 44cm Short Long	Under 44cm Over 44cm Short Long	Under 44cm Over 44cm Short Long	Under 44cm Over 44cm Short Long	Under 44cm Over 44cm Short Long	Under 44cm Over 44cm Short Long	Under 44cm Over 44cm Short Long

	Left	Right	Only Use c^{A1} if ordering a detached Juxta-Lite™ Ankle-Foot Wrap	
c^{A1} Arch Circ. (Foot Size)			19-26cm Small	24-31cm Med

Product Selection	
Legging Only	<input type="checkbox"/>
Legging with detached Compression Anklet ¹	<input type="checkbox"/>
Legging with detached Juxta-Lite™ Ankle-Foot Wrap	<input type="checkbox"/>

¹Note: The Juxta-Lite with 8-15mmHg Anklet will be available Spring, 2010.

Foot Piece (Ankle-Foot Wrap)

The Juxta-Lite™ legging can be ordered by itself, with a detached compression anklet, or with a detached Juxta-Lite Ankle-Foot Wrap. The Juxta-Lite Ankle-Foot Wrap can accommodate the circumference ranges listed above (**c^{A1}**), and a heel-to-base-of-toe length of 16cm or greater. If foot measurements fall outside of these ranges and patient fits into a standard size legging please contact us for additional assistance with your order.